

4
Linear
actuators



DC-04
Blood Donor Chair

Colour
Options



Actuators With Min IP X4
Protection For Quiet
Operation

Hand Control Unit



Control Units (Standard)
Battery (Optional)

*Control units pictures are
representative.*



70° Backrest Movement



80° Legrest Movement



57,5-80cm Height Movement



30 DegreeTrend Movement

Width (Overall)	Lenght (Overall)	Height Min.	Height Max.	Trend Angle	Sitting Surface	Backrest Angle	Legrest Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
76,5 cm	171,5 cm	57,5 cm	80 cm	30°	55x48,5 cm	70°	80°	230 kg	220 – 240 V AC	4

Arm supports working
parallel to the floor

Three-piece
upholstery with
poliurethan filling

ABS plastic
bottom chassis cover

Nano ceramic coated
electrostatic powder
painted metal parts

Bottom chassis profiles:
30x30x2mm
30x50x2mm

Mechanism profiles:
30x40x2mm

Upper chassis
profiles:
20x20x2mm

Material Features



B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr