



B K - 01
Baby Cot

GENERAL FEATURES



Width	Lenght	Height Min.	Height Max.	Mattress	Trendelenburg Angle
65 cm	45 cm	87 cm	40 cm	30x60x5 cm 28DNS	14°



B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.No:19/A1
Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr