




2
Lifting
columns

+

2
Linear
actuators



M2-DXS
Delivery Bed





Actuators With Min
IPX4 Protection
For Quiet Operation

Ø150 mm 360 ° Swivel
Central Brake Castors



Nurse Panel With
Function Locking Feature



Control Box
Built-in Battery



Squatting Bar



Waste Basin



Height Adjustable Iv Pole



Calf Support




Control units pictures are representative.




Siderails Control Panel

Sticker options
in different colours
for Head/Footboards
and side rails



X-Ray Cassete Holder


Underbed Light



Hand controller

Foot control







64-103cm Height Movement



12° Trendelenburg/Antitrendelenburg




70° Backrest Movement



Leg Extension

Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Lying Surface Closed	Lying Surface Open	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
103,5 cm	209,5 cm	64 cm	103 cm	85x144 cm	85x204 cm	85x142x12 cm 45x77x16 cm 28DNS	70°	30°	12°	12°	230 kg	220 – 240 V AC	4





B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr