


4
Linear
actuators



PLUS A8

Electromechanic Hospital Bed


Colour
Options





Actuators With Min
IPX4 Protection
For Quiet Operation

Ø150 mm 360 ° Swivel
Central Brake Castors



Hand Control Unit
Nurse Panel With
Function Locking Feature



Control Box
Built-in Battery




Manual CPR Arm

Mechanical Angle Indicator
on Siderails (backrest- trend angle)



Height Adjustable Iv Pole




ABS Plastic Bottom Chassis

Control units pictures are representative.



Siderails Control Panel


Underbed Light



Sticker options
in different colours
for Head/Footboards
and side rails



HPL Lying Surface



X-Ray Cassete Holder


Leg Extension Part




File Holder




Double Wheel




40-73cm Height Movement



12° Trendelenburg/Antitrendelenburg




70° Backrest Movement



30° Legrest Movement

Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Lying Surface	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
103,5 cm	212,5 cm	40 cm	73 cm	85x190 cm	85x190x12 cm 28DNS	70°	30°	12°	12°	230 kg	220 – 240 V AC	4





B. Kayacık Osb. Mahallesi
T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.
No:19/A1 Selçuklu / Konya / Turkey
Tel : +90 332 239 10 80
Fax : +90 332 239 10 82 (pbx)
E-Mail : bilgi@medikal2000.com.tr
www.medikal2000.com.tr