



**T-700
HYDRAULIC
EMERGENCY
STRETCHER**

- Nano ceramic coated electrostatic powder coated metal parts
- 4-part lying surface made of compact material
- Easy to open and close side rails
- Crash bumpers on all four corners
- Telescopic adjustable IV pole
- Plastic coated foldable push bars
- 360° swivel central locking wheels (200 mm)
- 5.Th Castor (Provides ease of movement in the desired direction on the stretcher)
- Hydraulically controlled height, trendelenburg and reverse trendelenburg adjustment with foot pedals
- Gas spring controlled back and foot adjustment
- Ability to adjust knee movements with gradual mechanism
- Aluminum telescopic columns
- ABS plastic coated sub-chassis
- Oxygen cylinder holder
- Patient restraint belt
- 28 density mattress (66x193x8cm)
- X-Ray Cassette carrier




**T-700
HYDRAULIC
EMERGENCY
STRETCHER**



1

Ayarlanabilir serum askısı
Height adjustable IV pole



2

Ayak Pedalı ile kolay Kullanım
Foot Pedal For Easy Use



3

5. Teker
5.Th Castor



4

ABS Alt Şase ve Merkezi Kilit Sistemi
ABS Plastic Bottom Chassis and Central Brake System



5

Katlanabilir İtme Barı
Foldable Push Bar



6

Sırt Kısımında Çift Amortisör
Sırt Kısımında Çift Amortisör



**T-700
HYDRAULIC
EMERGENCY
STRETCHER**



**FOOTREST
POSITION**



**BACKREST
POSITION**



**HEIGHT ADJUST.
POSITION**



**TRENDELENBURG
POSITION**



**R. TRENDELENBURG
POSITION**

**T-700
HYDRAULIC
EMERGENCY
STRETCHER**



General Features

Width (Over all)	Lenght (Over all)	Height Min	Height Max	Lying Surface	Backrest Angle	Legrest Angle	Trendelenburg Angle	R.trendelenburg Angle	Safety Working Capacity
84 cm	208,5 cm	58 cm	86 cm	60x193 cm	90	30	10	10	230 kgs



Medikal 2000

Technology for Health

B. Kavacık Osb. Mahallesi

T. Ziyaeddin Akbulut Cad. Fabrika Üretim Apt.No:19/A1

Selçuklu / Konya / Turkey

Tel : +90 332 239 10 80

Fax : +90 332 239 10 82 (pbx)

E-Mail : bilgi@medikal2000.com.tr

www.medikal2000.com.tr